

Pelvic Floor / Beckenboden / Plancher pelvien / Suelo Pélvico / Pavimento pelvico

FR	Douleur pelvien		DE	Becken. Schmerz				
IT	Dolore del pavimento pelvico		SP	Dolor en el suelo pélvico				
EN	Pelvic Floor Pain		▼1	▼2	▼3	▼4	▼5	
URO 1		20:00	STIM CON					
		[mm:ss]	20:00					
		[Hz]	3					
		[µS]	150					
		[sec]	1					
		[sec]	0					
		[sec]						
		[sec]						
		[sec]						

FR	Incontinence par urgence 1		DE	Dranginkontinenz 1				
IT	Incontinenza urgenza 1		SP	Incontinencia imperiosa 1				
EN	Urge Incontinence 1		▼1	▼2	▼3	▼4	▼5	
URO 2		20:00	STIM SYN					
		[mm:ss]	20:00					
		[Hz]	10					
		[µS]	250					
		[sec]	1					
		[sec]	0					
		[sec]	5					
		[sec]	5					
		[sec]						

FR	Incontinence d'effort 1		DE	Stressinkontinenz 1				
IT	Incontinenza urinaria da stress 1		SP	Incontinencia de esfuerzo 1				
EN	Stress Incontinence 1		▼1	▼2	▼3	▼4	▼5	
URO 3		20:00	STIM SYN					
		[mm:ss]	20:00					
		[Hz]	40					
		[µS]	200					
		[sec]	1					
		[sec]	0					
		[sec]	6					
		[sec]	15					
		[sec]						

FR	Incontinence d'effort 2		DE	Stressinkontinenz 2				
IT	Incontinenza urinaria da stress 2		SP	Incontinencia de esfuerzo 2				
EN	Stress Incontinence 2		▼1	▼2	▼3	▼4	▼5	
URO 4		20:00	STIM SYN					
		[mm:ss]	20:00					
		[Hz]	30					
		[µS]	200					
		[sec]	1					
		[sec]	0					
		[sec]	5					
		[sec]	8					
		[sec]						

FR	Fréquence/Urgence 1		DE	Häufigkeit /Drang 1				
IT	Frequenza/Stimolo 1		SP	Frecuencia/Imperiosa 1				
EN	Frequency/Urge 1		▼1	▼2	▼3	▼4	▼5	
URO 5		20:00	STIM SYN					
		[mm:ss]	20:00					
		[Hz]	10					
		[µS]	200					
		[sec]	1					
		[sec]	0					
		[sec]	5					
		[sec]	5					
		[sec]						

FR	Fréquence/Urgence 2		DE	Häufigkeit /Drang 2				
IT	Frequenza/Stimolo 2		SP	Frecuencia/Imperiosa 2				
EN	Frequency/Urge 2	▼1	▼2	▼3	▼4	▼5		
URO 6		🕒 15:00	STIM CON					
		[mm:ss]	15:00					
		[Hz]	10					
		[µS]	200					
		[sec]						
		[sec]						
		[sec]						
		[sec]						
		[sec]						

FR	Fréquence/Urgence 3		DE	Häufigkeit /Drang 3				
IT	Frequenza/Stimolo 3		SP	Frecuencia/Imperiosa 3				
EN	Frequency/Urge 3	▼1	▼2	▼3	▼4	▼5		
URO 7		🕒 20:00	STIM CON					
		[mm:ss]	20:00					
		[Hz]	10					
		[µS]	220					
		[sec]						
		[sec]						
		[sec]						
		[sec]						
		[sec]						

FR	Manque de sensation		DE	Sensibilitätsverlust				
IT	Mancanza di sensibilità		SP	Ausencia de sensibilidad				
EN	Lack of Sensitivity	▼1	▼2	▼3	▼4	▼5		
URO 8		🕒 25:00	STIM SYN	STIM SYN	STIM SYN	STIM SYN	STIM SYN	
		[mm:ss]	03:00	10:00	05:00	04:00	03:00	
		[Hz]	3	10	20	30	40	
		[µS]	250	250	250	200	200	
		[sec]	0.8	0.8	0.8	0.7	0.7	
		[sec]	0	0	0	0	0	
		[sec]	4	4	4	4	4	
		[sec]	4	4	4	6	6	
		[sec]						

NeuroTrac® MyoPlus 2/4 Programmes

FR	Exercice – Prophylaxie d’une légère faiblesse musculaire		DE	Work Out - Prophylaxe von leichter Muskelschwäche			
IT	Esercizio - profilassi di lieve debolezza muscolare		SP	Ejercicio - Profilaxis de la debilidad muscular leve			
EN	Work Out - Prophylaxis of mild muscular weakness		▼1	▼2	▼3	▼4	▼5
URO 9		⌚ 45:00	STIM SYN	STIM SYN	STIM SYN	STIM SYN	STIM SYN
		[mm:ss]	04:00	15:00	08:00	08:00	10:00
		[Hz]	20	10	20	30	10
		[µS]	250	250	250	200	250
		[sec]	0.8	0.8	0.8	0.6	0.8
		[sec]	0	0	0	0	0
		[sec]	5	5	5	5	5
		[sec]	5	5	5	6	7
		[sec]					
FR	Accroissement de l’endurance		DE	Aufbau von Ausdauer			
IT	Sviluppo della resistenza a lungo termine alla fatica		SP	Aumentar el aguante y la resistencia ante la fatiga			
EN	Building Up Endurance		▼1	▼2	▼3	▼4	▼5
URO 10		⌚ 20:00	STIM SYN				
		[mm:ss]	20:00				
		[Hz]	20				
		[µS]	250				
		[sec]	0.8				
		[sec]	0				
		[sec]	5				
		[sec]	5				
		[sec]					
FR	Relâchement du muscle pelvien		DE	Entspannung des Beckenmuskels			
IT	Rilassamento la muscolatura pelvica		SP	Como relajar el músculo pélvico			
EN	Relaxing the Pelvic M.		▼1	▼2	▼3	▼4	▼5
URO 11		⌚ 20:00	STIM SYN				
		[mm:ss]	20:00				
		[Hz]	2				
		[µS]	220				
		[sec]	1.2				
		[sec]	1.2				
		[sec]	6				
		[sec]	10				
		[sec]					

FR	Contenance après l'accouchement			DE	Kontinenz nach der Geburt		
IT	Continenza dopo il parto			SP	Incont. después de un parto		
EN	Contenance after Childbirth		▼1	▼2	▼3	▼4	▼5
URO 12		🕒 20:00	EMG	ETS	EMG	STIM SYN	
		[mm:ss]	02:10	04:15	01:00	05:00	
		[Hz]		40		30	
		[µS]		220		220	
		[sec]		1		1	
		[sec]		0.8		0.8	
		[sec]	5	5/6	4	5	
		[sec]	8	12	8	12	
		[sec]	10	15	5		

FR	Incontinence d'effort 3			DE	Stress Incontinenz 3		
IT	Incont. urinaria da stress 3			SP	Incontinencia de esfuerzo 3		
EN	Stress Incontinence 3		▼1	▼2	▼3	▼4	▼5
URO 13		🕒 15:49	EMG	STIM SYN	EMG	STIM SYN	EMG
		[mm:ss]	01:05	08:00	00:48	05:00	00:56
		[Hz]		35		30	
		[µS]		220		220	
		[sec]		0.8		1	
		[sec]		0.8		0.8	
		[sec]	5	5	4	5	4
		[sec]	8	7	8	12	10
		[sec]	5		4		4

FR	Incontinence d'effort 4			DE	Stress Incontinenz 4		
IT	Incont. urinaria da stress 4			SP	Incontinencia de esfuerzo 4		
EN	Stress Incontinence 4		▼1	▼2	▼3	▼4	▼5
URO 14		🕒 11:39	EMG	STIM SYN	EMG	ETS	EMG
		[mm:ss]	01:05	08:00	00:48	00:50	00:56
		[Hz]		40		35	
		[µS]		220		220	
		[sec]		0.8		0.8	
		[sec]		0.8		0.8	
		[sec]	5	5	4	5/5	4
		[sec]	8	7	8	5	10
		[sec]	5		4	5	4

FR	Incontinence par urgence 2			DE	Dranginkontinenz 2			
IT	Incontinenza urgenza 2			SP	Incontinenca imperiosa 2			
EN	Urge Incontinence 2			▼1	▼2	▼3	▼4	▼5
URO 15		11:36	EMG	STIM SYN	EMG	STIM SYN	EMG	
		[mm:ss]	00:48	05:00	00:48	04:00	01:00	
		[Hz]		10		10		
		[µS]		220		220		
		[sec]		0.8		0.8		
		[sec]		0.8		0.8		
		[sec]	4	5	4	4	5	
		[sec]	8	7	8	7	10	
		[sec]	4		4		4	
FR	Incontinence par urgence 3			DE	Dranginkontinenz 3			
IT	Incontinenza urgenza 3			SP	Incontinenca imperiosa 3			
EN	Urge Incontinence 3			▼1	▼2	▼3	▼4	▼5
URO 16		11:36	EMG	STIM SYN	EMG	STIM SYN	EMG	
		[mm:ss]	00:48	05:00	00:48	04:00	01:00	
		[Hz]		10		10		
		[µS]		220		220		
		[sec]		0.8		0.8		
		[sec]		0.8		0.8		
		[sec]	4	7	4	5	5	
		[sec]	8	5	8	5	10	
		[sec]	4		4		4	
FR	Incontinence par urgence 4 - maintien			DE	Dranginkontinenz 4 - Aufrechterhaltung			
IT	Incontinenza urgenza 4 - mantenimento			SP	Incontinenca imperiosa 4 - Mantenimiento			
EN	Urge Incontinence 4 - maintenance			▼1	▼2	▼3	▼4	▼5
URO 17		17:00	EMG	STIM SYN	EMG			
		[mm:ss]	01:12	15:00	00:48			
		[Hz]		10				
		[µS]		240				
		[sec]		0.8				
		[sec]		0.8				
		[sec]	4	5	4			
		[sec]	8	7	8			
		[sec]	6		4			

FR	Muscles distendus du plancher pelvien	DE	Schlaffe Beckenbodenmuskulatur				
IT	Flaccidità dei muscoli del pavimento pelvico	SP	Flacidez muscular en el suelo pélvico				
EN	Flaccid Pelvic Floor Muscles	▼1	▼2	▼3	▼4	▼5	
URO 18		5:40	STIM SYN	ETS			
		[mm:ss]	05:00	00:40			
		[Hz]	10	20			
		[µS]	250	250			
		[sec]	0.8	0.8			
		[sec]	0.8	0.8			
		[sec]	6	5/6			
		[sec]	7	5			
		[sec]		4			

FR	Faiblesse des muscles du plancher pelvien	DE	Schwache Beckenbodenmuskulatur				
IT	Debolezza dei muscoli del pavimento pelvico	SP	Debilidad muscular en el suelo pélvico				
EN	Weak Pelvic Floor Muscles	▼1	▼2	▼3	▼4	▼5	
URO 19		7:12	EMG	ETS	EMG	ETS	STIM SYN
		[mm:ss]	00:48	00:44	00:56	00:44	04:00
		[Hz]		10		20	30
		[µS]		220		220	220
		[sec]		0.8		0.7	0.8
		[sec]		0.8		0.8	0.8
		[sec]	5	5/6	4	4/6	4
		[sec]	7	6	10	7	7
		[sec]	4	4	4	4	

SPORTS / DEPORTE










FR	Echauffement		DE	Aufwärmen				
IT	Riscaldamento		SP	Calentamiento				
EN	Warm Up		▼1	▼2	▼3	▼4	▼5	
SPORTS 1		30:00	STIM SYN	STIM MOD				
		[mm:ss]	20:00	10:00				
		[Hz]	10	15 - 2				
		[µS]	220	220				
		[sec]	0.8	0.8				
		[sec]	0.4	0.4				
		[sec]	5	5				
		[sec]	6	6				
		[sec]						










FR	Récupération active		DE	Aktive Erholung				
IT	Recupero attivo		SP	Recuperación activa				
EN	Active Recovery		▼1	▼2	▼3	▼4	▼5	
SPORTS 2		39:00	STIM SYN	STIM SYN	STIM SYN	STIM SYN		
		[mm:ss]	15:00	10:00	07:00	07:00		
		[Hz]	10	8	5	2		
		[µS]	250	250	220	220		
		[sec]	0.8	0.8	0.8	0.8		
		[sec]	0.4	0.4	0.4	0.4		
		[sec]	5	5	5	5		
		[sec]	7	6	6	5		
		[sec]						










FR	Endurance		DE	Ausdauertraining				
IT	Durata		SP	Aguante				
EN	Endurance		▼1	▼2	▼3	▼4	▼5	
SPORTS 3		45:00	STIM CON	STIM MOD	STIM MOD	STIM MOD		
		[mm:ss]	05:00	15:00	20:00	05:00		
		[Hz]	5	15 - 5	20 - 5	15 - 2		
		[µS]	250	250	220	175 - 250		
		[sec]		2	2	4		
		[sec]		2	2	4		
		[sec]		5	6	8		
		[sec]		8	8	8		
		[sec]						

FR	Résistance 1		DE	Ausdauer 1				
IT	Resistenza 1		SP	Resistencia 1				
EN	Resistance 1		▼1	▼2	▼3	▼4	▼5	
SPORTS 4		35:00	STIM CON	STIM MOD	STIM SYN			
		[mm:ss]	05:00	20:00	10:00			
		[Hz]	5	20 - 2	10			
		[µS]	250	240	250			
		[sec]		3	0.8			
		[sec]		2	0.4			
		[sec]		5	5			
		[sec]		10	5			
		[sec]						

FR	Résistance 2		DE	Ausdauer 2				
IT	Resistenza 2		SP	Resistencia 2				
EN	Resistance 2		▼1	▼2	▼3	▼4	▼5	
SPORTS 5		35:00	STIM CON	STIM MOD	STIM SYN			
		[mm:ss]	05:00	20:00	10:00			
		[Hz]	5	40 - 2	10			
		[µS]	250	240	250			
		[sec]		3	0.8			
		[sec]		2	0.4			
		[sec]		5	5			
		[sec]		10	5			
		[sec]						

FR	Résistance 3		DE	Ausdauer 3				
IT	Resistenza 3		SP	Resistencia 3				
EN	Resistance 3			▼1	▼2	▼3	▼4	▼5
SPORTS 6  35:00			STIM CON	STIM MOD	STIM SYN			
		[mm:ss]	05:00	20:00	10:00			
		[Hz]	5	50 - 2	10			
		[µS]	250	240	250			
		[sec]		3	0.8			
		[sec]		2	0.4			
		[sec]		5	5			
		[sec]		10	6			
		[sec]						

FR	Force et résistance 1		DE	Kraft/Ausdauer 1				
IT	Resistenza Forza 1		SP	Resistencia fuerza 1				
EN	Resistance Strength 1			▼1	▼2	▼3	▼4	▼5
SPORTS 7  35:00			STIM CON	STIM MOD	STIM SYN			
		[mm:ss]	05:00	20:00	10:00			
		[Hz]	5	60 - 2	10			
		[µS]	250	250	240			
		[sec]		3	0.8			
		[sec]		2	0.4			
		[sec]		5	5			
		[sec]		10	7			
		[sec]						

FR	Force et résistance 2		DE	Kraft/Ausdauer 2				
IT	Resistenza Forza 2		SP	Resistencia fuerza 2				
EN	Resistance Strength 2			▼1	▼2	▼3	▼4	▼5
SPORTS 8  35:00			STIM CON	STIM MOD	STIM SYN			
		[mm:ss]	05:00	20:00	10:00			
		[Hz]	5	75 - 2	10			
		[µS]	250	250	240			
		[sec]		3	0.8			
		[sec]		2	0.6			
		[sec]		5	5			
		[sec]		10	8			
		[sec]						

FR	Force et résistance 3		DE	Kraft/Ausdauer 3				
IT	Resistenza Forza 3		SP	Resistencia fuerza 3				
EN	Resistance Strength 3		▼1	▼2	▼3	▼4	▼5	
SPORTS 9		35:00	STIM CON	STIM MOD	STIM SYN			
		[mm:ss]	05:00	20:00	10:00			
		[Hz]	5	90 - 2	10			
		[µS]	250	250	240			
		[sec]		3	0.8			
		[sec]		2	0.6			
		[sec]		5	5			
		[sec]		10	8			
		[sec]						

FR	Force 1		DE	Kraft 1				
IT	Forza 1		SP	Fuerza 1				
EN	Strength 1		▼1	▼2	▼3	▼4	▼5	
SPORTS 10		35:00	STIM CON	STIM MOD	STIM SYN			
		[mm:ss]	05:00	20:00	10:00			
		[Hz]	5	40 - 2	10			
		[µS]	250	290	240			
		[sec]		2	2			
		[sec]		2	2			
		[sec]		6	5			
		[sec]		20	7			
		[sec]						

FR	Force 2		DE	Kraft 2				
IT	Forza 2		SP	Fuerza 2				
EN	Strength 2		▼1	▼2	▼3	▼4	▼5	
SPORTS 11		35:00	STIM CON	STIM MOD	STIM MOD			
		[mm:ss]	05:00	20:00	10:00			
		[Hz]	5	50 - 5	10 - 5			
		[µS]	250	340	240			
		[sec]		2	2			
		[sec]		2	2			
		[sec]		6	5			
		[sec]		20	7			
		[sec]						

FR	Force 3		DE	Kraft 3				
IT	Forza 3		SP	Fuerza 3				
EN	Strength 3		▼1	▼2	▼3	▼4	▼5	
SPORTS 12		35:00	STIM CON	STIM MOD	STIM MOD			
		[mm:ss]	05:00	20:00	10:00			
		[Hz]	5	60 - 5	10 - 5			
		[µS]	250	375	240			
		[sec]		2	2			
		[sec]		2	2			
		[sec]		6	5			
		[sec]		22	8			
		[sec]						









FR	Force max. 1		DE	Max. Kraft 1				
IT	Forza Massima 1		SP	Fuerza Máx. 1				
EN	Max Strength 1		▼1	▼2	▼3	▼4	▼5	
SPORTS 13		35:00	STIM CON	STIM MOD	STIM MOD			
		[mm:ss]	05:00	20:00	10:00			
		[Hz]	5	70 - 5	10 - 5			
		[µS]	250	390	240			
		[sec]		2	2			
		[sec]		2	2			
		[sec]		6	5			
		[sec]		22	8			
		[sec]						









FR	Force max. 2		DE	Max. Kraft 2				
IT	Forza Massima 2		SP	Fuerza Máx. 2				
EN	Max Strength 2		▼1	▼2	▼3	▼4	▼5	
SPORTS 14		35:00	STIM CON	STIM MOD	STIM MOD			
		[mm:ss]	05:00	20:00	10:00			
		[Hz]	5	80 - 5	10 - 5			
		[µS]	250	390	240			
		[sec]		2	2			
		[sec]		2	2			
		[sec]		6	5			
		[sec]		22	8			
		[sec]						









FR	Force max. 3		DE	Max. Kraft 3				
IT	Forza Massima 3		SP	Fuerza Máx. 3				
EN	Max Strength 3		▼1	▼2	▼3	▼4	▼5	
SPORTS 15		35:00	STIM CON	STIM MOD	STIM MOD			
		[mm:ss]	05:00	20:00	10:00			
		[Hz]	5	90 - 5	10 - 5			
		[µS]	250	400	240			
		[sec]		2	2			
		[sec]		2	2			
		[sec]		6	5			
		[sec]		25	8			
		[sec]						










FR	Force max. 4		DE	Max. Kraft 4				
IT	Forza Massima 4		SP	Fuerza Máx. 4				
EN	Max Strength 4		▼1	▼2	▼3	▼4	▼5	
SPORTS 16		35:00	STIM CON	STIM MOD	STIM MOD			
		[mm:ss]	05:00	20:00	10:00			
		[Hz]	5	90 - 5	10 - 5			
		[µS]	250	400	240			
		[sec]		2	2			
		[sec]		2	2			
		[sec]		6	5			
		[sec]		25	8			
		[sec]						

FR	Force explosive 1		DE	Schnellkraft 1				
IT	Forza esplosiva 1		SP	Fuerza explosiva 1				
EN	Explosive Force 1		▼1	▼2	▼3	▼4	▼5	
SPORTS 17		35:00	STIM CON	STIM MOD	STIM MOD			
		[mm:ss]	05:00	20:00	10:00			
		[Hz]	10	60 - 5	10 - 5			
		[µS]	250	200 - 400	240			
		[sec]		2	2			
		[sec]		2	2			
		[sec]		7	5			
		[sec]		15	10			
		[sec]						

FR	Force explosive 2		DE	Schnellkraft 2				
IT	Forza esplosiva 2		SP	Fuerza explosiva 2				
EN	Explosive Force 2		▼1	▼2	▼3	▼4	▼5	
SPORTS 18		🕒 35:00	STIM CON	STIM MOD	STIM MOD			
		[mm:ss]	05:00	20:00	10:00			
		[Hz]	10	70 - 5	10 - 5			
		[µS]	250	200 - 400	240			
		[sec]		2	2			
		[sec]		2	2			
		[sec]		7	5			
		[sec]		15	10			
		[sec]						

FR	Résistance à la fatigue musculaire		DE	Ermüdung				
IT	Resistenza alla fatica muscolare		SP	Resistencia a la fatiga muscular				
EN	Fatigue resistance		▼1	▼2	▼3	▼4	▼5	
SPORTS 19		🕒 30:00	STIM CON	STIM MOD	STIM CON			
		[mm:ss]	05:00	15:00	10:00			
		[Hz]	10	40 - 5	5			
		[µS]	240	275	220			
		[sec]		3				
		[sec]		2				
		[sec]		6				
		[sec]		20				
		[sec]						

FR	Antagonistes		DE	Agonisten Stärkung				
IT	Rafforzamento dei muscoli antagonisti		SP	Fortalecimiento de músculos agonistas				
EN	Agonistic strengthening		▼1	▼2	▼3	▼4	▼5	
SPORTS 20		🕒 10:00	STIM CON	STIM MOD				
		[mm:ss]	06:00	04:00				
		[Hz]	5	30 - 2				
		[µS]	220	220				
		[sec]		1				
		[sec]		0.6				
		[sec]		8				
		[sec]		20				
		[sec]						

FR	Soulagement de la douleur		DE	Schmerzlinderung			
IT	Antidolorifico		SP	Alivio del dolor			
EN	Pain relief		▼1	▼2	▼3	▼4	▼5
SPORTS 21		🕒 14:00	STIM CON	STIM MOD	STIM CON		
		[mm:ss]	05:00	04:00	05:00		
		[Hz]	5	30 - 20	10		
		[µS]	220	250	220		
		[sec]		3			
		[sec]		3			
		[sec]		8			
		[sec]		20			
		[sec]					

REHABILITACION / REABILITAZIONE

FR	Contrôle des muscles posturaux		DE	Haltemuskulator Kontrolle				
IT	Controllo dei muscoli posturali		SP	Control de postura				
EN	Postural muscles control		▼1	▼2	▼3	▼4	▼5	
REHAB 1		32:00	STIM CON	STIM SYN	STIM SYN			
		[mm:ss]	02:00	25:00	05:00			
		[Hz]	5	10	10			
		[µS]	220	250	220			
		[sec]		2	1			
		[sec]		0.6	0.6			
		[sec]		8	5			
		[sec]		10	7			
		[sec]						

FR	Colonne vertébrale		DE	Wirbelsäule				
IT	Colonna vertebrale		SP	Columna vertebral				
EN	Spinal Column		▼1	▼2	▼3	▼4	▼5	
REHAB 2		30:00	STIM CON					
		[mm:ss]	30:00					
		[Hz]	70					
		[µS]	200					
		[sec]						
		[sec]						
		[sec]						
		[sec]						
		[sec]						

FR	Parésie des extrémités		DE	Paretische Extremitäten				
IT	Paresi delle estremità		SP	Parálisis de las extremidades				
EN	Paretic Extremities		▼ 1	▼ 2	▼ 3	▼ 4	▼ 5	
REHAB 3		15:00	STIM SYN					
		[mm:ss]	15:00					
		[Hz]	20					
		[µS]	250					
		[sec]	4					
		[sec]	1.5					
		[sec]	10					
		[sec]	30					
		[sec]						

FR	Paralysie faciale 1		DE	Fazialisparese 1				
IT	Paralisi facciale 1		SP	Parálisis facial 1				
EN	Facial Palsy 1		▼ 1	▼ 2	▼ 3	▼ 4	▼ 5	
REHAB 4		40:00	STIM SYN					
		[mm:ss]	40:00					
		[Hz]	8					
		[µS]	80					
		[sec]	1					
		[sec]	0.6					
		[sec]	6					
		[sec]	10					
		[sec]						

FR	Paralysie faciale 2		DE	Fazialisparese 2				
IT	Paralisi facciale 2		SP	Parálisis facial 2				
EN	Facial Palsy 2		▼ 1	▼ 2	▼ 3	▼ 4	▼ 5	
REHAB 5		40:00	STIM SYN					
		[mm:ss]	40:00					
		[Hz]	20					
		[µS]	80					
		[sec]	1					
		[sec]	0.6					
		[sec]	6					
		[sec]	10					
		[sec]						

FR	Atrophie		DE	Muskelatrophie				
IT	Atrofia da disuso		SP	Atrofia por inactividad				
EN	Disuse Atrophy	▼1	▼2	▼3	▼4	▼5		
REHAB 6		60:00	STIM SYN					
		[mm:ss]	60:00					
		[Hz]	10					
		[µS]	220					
		[sec]	1					
		[sec]	0.6					
		[sec]	6					
		[sec]	7					
		[sec]						

FR	Programme 1 suite à un AVC		DE	Nach dem Schlaganfall 1				
IT	Dopo ictus 2		SP	Después de un ictus 1				
EN	After Stroke 1	▼1	▼2	▼3	▼4	▼5		
REHAB 7		60:00	STIM SYN					
		[mm:ss]	60:00					
		[Hz]	10					
		[µS]	220					
		[sec]	2					
		[sec]	1					
		[sec]	6					
		[sec]	10					
		[sec]						

FR	Programme 2 suite à un AVC		DE	Nach dem Schlaganfall 1				
IT	Dopo ictus 2		SP	Después de un ictus 2				
EN	After Stroke 2		▼1	▼2	▼3	▼4	▼5	
REHAB 8	60:00	STIM SYN	STIM SYN					
	[mm:ss]	50:00	10:00					
	[Hz]	10	40					
	[µS]	220	220					
	[sec]	2	2					
	[sec]	1	1					
	[sec]	6	6					
	[sec]	8	8					
	[sec]							

FR	Sclérose en plaques		DE	Multiple Sklerose				
IT	Sclerosi a placche		SP	Esclerosis múltiple				
EN	Multiple Sclerosis		▼1	▼2	▼3	▼4	▼5	
REHAB 9	25:00	STIM SYN	STIM SYN					
	[mm:ss]	10:00	15:00					
	[Hz]	3	10					
	[µS]	220	220					
	[sec]	2	2					
	[sec]	1	1					
	[sec]	6	6					
	[sec]	8	10					
	[sec]							